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ULTIMA REVISIONE: 01/2017

The Wave Innovation Group srls
Via Barucchi, 47/A - 37139 Verona (Italy)
www.thewaveresearch.it - info@thewaveresearch.it
www.thewaveinnovation.com - info@thewaveinnovation.com

Proage

HA Biostimulating Reparative Phase

DESCRIPTION:

PROAGE is a sterile, biodegradable, isotonic, intradermal gel. PROAGE consists of medium chain ($1.0\text{--}1.5 \times 10^6$ Dalton) hyaluronic acid, obtained from *Streptococcus equi* bacteria, formulated to a concentration of 20 mg/ml in a physiologic buffer. Each box contains one syringe of PROAGE, two 30G ½" disposable sterile needles reserved for injection of PROAGE and a product leaflet. A set of two labels showing the batch number is contained in the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

COMPOSITION:

Sodium hyaluronate 20 mg/ml, sodium chloride, sodium dihydrogen phosphate dihydrate, dibasic sodium phosphate dodecahydrate, WFI grade water. One syringe contains 2.0 ml of non-pyrogenic gel, sterilised using moist heat.

INDICATIONS:

PROAGE is indicated:

- in the physiological process of ageing of the skin, the effects of which include insipidation of the horny layer (roughness and reduced radiance of the skin) and changes in the elastic fibres of the dermis (wrinkles),
- in the dermal tissue repair process, in cases of scar results following superficial cutaneous trauma (e.g. acne and chickenpox scars).

EXCLUSION CRITERIA:

PROAGE must not be used in:

- patients who tend to develop hypertrophic scarring,
- patients with a history of autoimmune disease or who are receiving immune therapy,
- patients who are known to be hypersensitive to hyaluronic acid,
- pregnant or breastfeeding women,
- patients under 18 years of age.

Anticoagulated patients or patients receiving platelet aggregation inhibitors should not be treated with PROAGE without consulting their doctors.

PROAGE must not be used in areas presenting cutaneous, inflammatory and/or infectious processes (e.g. acne, herpes ...). PROAGE must not be used in association with laser therapy, chemical peeling or dermabrasion.

PRECAUTIONS FOR USE:

PROAGE is only indicated for intradermal injections. Sensitive skin may be pretreated using a local anaesthetic patch or cream. Please note that any anaesthesia may cause redness or local hypersensitivity. There are no available clinical data (efficiency, tolerance) about injecting PROAGE into an area which has already been treated with another filling product. Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV or using saunas or Turkish baths for one week after the injection. If the needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle.

Do not inject into blood vessels. Do not use PROAGE in bones, tendons, ligaments or muscles. Do not inject PROAGE into naevi. Discard the syringe and remaining product after use.

Do not use if packaging has been damaged.

INTERACTIONS WITH OTHER AGENTS:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore PROAGE should never be placed in contact with these substances or with medical-surgical instruments that have been in contact with these substances.

UNDESIRE SIDE EFFECTS:

Physicians must inform the patient that there are potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list):

- Inflammatory reactions (e.g. redness, oedema, erythema ...) which may be associated with itching and pain on pressure after the injection. These reactions may last for a week.
- Indurations or nodules at the injection site especially in cases of too superficial placement.
- Discolouration of the injection site.
- Poor effect if PROAGE is injected incorrectly.

Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other secondary effect which develops. The physician should treat these as appropriate. Any other undesirable side effects associated with the injection of PROAGE must be reported to the distributor and/or to the manufacturer.

METHODS OF USE:

PROAGE must be injected into non-inflamed, disinfected, healthy skin. The technique used is essential for the success of the treatment. Therefore this device must only be used by doctors who have received specific training on the injection technique. The area to be treated should be disinfected thoroughly prior to the injection. Use the 30G ½" needle which is provided with the syringe and inject slowly by applying the appropriate injection technique.

Inject PROAGE at room temperature and with strict asepsis conditions.

Inject PROAGE intradermally with a linear technique or with picotage at a medium deep level.

The amount injected will depend on the wrinkle or depression which is to be corrected. After the injection, doctors may apply a light massage in order to distribute the product uniformly.


An initial course of three treatment sessions at one week intervals is recommended, followed if necessary by monthly maintenance sessions.

WARNINGS:


Verify the integrity of the syringe and the expiry date before use. Do not use any other needle or syringe than provided by the manufacturer. Do not re-use; quality and sterility can only be guaranteed for an originally closed syringe. The re-use of the product creates a potential infection risk for patients or users. Do not re-sterilise. After use the needle must be disposed of in a suitable container.

INSTRUCTIONS FOR ASSEMBLY OF THE SYRINGE NEEDLE:

A. Carefully unscrew the cap of the tip of the syringe, being particularly careful to avoid contact with the opening.

B. Gently grip the needle guard and mount the needle on the luer-lock mount, screwing it tight until a slight counter-pressure is felt in order to ensure an airtight grip and prevent leakage of the liquid during administration.

STORAGE:

Store PROAGE at 2–25°C (36–77°F) in a dry place in the original box. Protect from light, heat and frost and handle with care.

SYMBOLS

| | | | |
|--|---------------------------------|--|----------------------------------|
| | Consult instructions for use | | Store between +2°C and +25°C |
| | Do not re-sterilise | | Do not use if package is damaged |
| | Sterilised using steam | | Do not reuse |
| | Sterilised using ethylene oxide | | Manufacturer |
| | Keep dry | | Batch number |
| | Keep away from sunlight | | Use by |

ONLY FOR MEDICAL PRESCRIPTION

LAST REVISED: 01/2017

The Wave Innovation Group srls
Via Barucchi, 47/A - 37139 Verona (Italy)
www.thewaveresearch.it - info@thewaveresearch.it
www.thewaveinnovation.com - info@thewaveinnovation.com

Proage

HA Biostimulating Reparative Phase

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